

Installed Version Information

Application Software updated 30-Jan-2026 4:04:20 AM
Release 5.7.1.4

BOS For MR Console Computer updated 20-Feb-2020 12:02:55 PM
BOSWin10_H for Z420_20 Host computer

#####

Hardware

#

#

- Intel® Xeon® Quad-Core E5-1603 Processor 2.80 GHz CPU

#

- 32 GB DDR3 SDRAM

#

- On-board Intel® 82579LM Gigabit ethernet controller

#

- Integrated SATA 6 GB/s Controller with 2 SATA disks 10K RPM 500GB

#

#

- PCIe2 X4 Intel Pro Gigabit CT1 desktop adapter

#

- PCIe3 X16 NVIDIA Quadro Q410 display adapter

#

- PCIe2 x8 Free

#

- PCIe2 X8 Intel Pro Gigabit CT2 desktop adapter

#

- PCIe3 X8 Axxon LF734KB RS232

#

- PCIe3 X16 NVIDIA Quadro Q410 display adapter

#

- PCI free

#

#####

(or)

#####

Hardware

#

#

- Intel® Xeon® Six-Core E5-1650 CPU

#

- 64 GB DDR3 SDRAM

#

- On-board Intel® 82579LM Gigabit ethernet controller

#

- Integrated SATA 6 GB/s Controller with 1 SSD & 1 SATA disks 10K RPM 500GB #

#

- PCIe2 X4 Intel Pro Gigabit CT1 desktop adapter

#

- PCIe3 X16 NVIDIA Quadro Q410 display adapter

#

- PCIe2 x8 DNA NIC CARD

#

- PCIe2 X8 Free

#

- PCIe3 X8 Free

#

- PCIe3 X16 NVIDIA Quadro Q410 display adapter

#

- PCI BULK DATA CARD Interface

#

#####

OS

windows 10 Enterprise 2016 LTSC

System type Z420

Build Number 01.31.20.22

Date 2020-Feb-19

Commercial Package Description updated 11-May-2021 8:57:31 AM

File_name : NL9O_781002_50700_41480.lic

Fingerprint : 41480

Site_name : MOR INSIDE LTD

City : ISRAEL

Country : IL

System_type : ING-15M

Release : 50700

This license file supports the R5.7.x based SW releases.

Lock_type : 20171

Encryption_type : C

Generation_date : 16-Mar-2021 08:55:28

#-----

#-----

-

ENABLED COMMERCIAL OPTIONS:

| # COM_OPT | DESCRIPTION | CATEGORY |
|-----------|-------------|----------|
|-----------|-------------|----------|

| | | |
|----------------------|--|--|
| EXP_DATEORDER_NUMBER | | |
|----------------------|--|--|

| | | |
|-----------|---------------|----------------|
| # NMRF109 | 3D BreastView | |
| Permanent | | 1303462/000020 |

| | | |
|-----------|--------------|----------------|
| # NMRF088 | 3D NERVEVIEW | |
| Permanent | | 1303462/000020 |

| | | |
|-----------|---------------|----------------|
| # NMRF108 | 3D PelvisView | |
| Permanent | | 1303462/000020 |

| | | |
|-----------|--------------|----------------|
| # NMRF080 | 3D SpineView | |
| Permanent | | 1303462/000020 |

| | | |
|-----------|----------------|--|
| # NMRZ683 | ASL SPECIALIST | |
|-----------|----------------|--|

| | | |
|--------------|------------------------------------|----------------|
| Permanent | | 1303462/000020 |
| # NMRF064 | AUTOVOICE | |
| Permanent | | 1303462/000020 |
| # NMRF082 | BLACK BLOOD IMAGING | |
| Permanent | | 1303462/000020 |
| # NMRF117 | CARDIAC EXPERT SPEC | |
| Permanent | | 1303462/000020 |
| # NMRF116 | CORONARY ACQUISITION | |
| Permanent | | 1303462/000020 |
| # NMRF206 | Compressed SENSE Pack | |
| Permanent | | 1303462/000020 |
| # | including: | |
| # | NMRF203 Compressed SENSE Essential | |
| # | NMRF202 Compressed SENSE MSK | |
| # | NMRF204 Compressed SENSE Body | |
| # | NMRF205 Compressed SENSE Cardiac | |
| # CONTRCARDS | Contrast Cards | |
| Permanent | | 1303462/000020 |
| # NMRF075 | FIBERTRAK SPECIALIST | |
| Permanent | | 1303462/000020 |
| # INGENIA | INGENIA CONFIGURATION | |
| Permanent | | 1303462/000020 |
| # NMRF129 | MDIXON XD FFE MULTISTATION | |
| Permanent | | 1303462/000020 |
| # NMRF104 | MDIXON XD FFE SPECIALIST | |
| Permanent | | 1303462/000020 |
| # NMRF093 | MDIXON XD TSE SPECIALIST | |
| Permanent | | 1303462/000020 |
| # NMRZ686 | MOBILE SYSTEM | |
| Permanent | | 1303462/000020 |
| # NMRB552 | MR PROCESS OPTION | |
| Permanent | | 1303462/000020 |
| # NMRF079 | MULTIVANEXD | |
| Permanent | | 1303462/000020 |
| # NMRF095 | O-MAR XD SPECIALIST | |

| | | |
|-----------|------------------------------|----------------|
| Permanent | | 1303462/000020 |
| # NMRZ233 | R5.7 based order | |
| Permanent | | 1303462/000020 |
| # NMRB525 | RECONEXCEL | |
| Permanent | | 1303462/000020 |
| # NMRF027 | SCANTOOLS PREMIUM | |
| Permanent | | 1303462/000020 |
| # | including: | |
| # | - Scantools pro + | |
| # | - K-T Blast | |
| # | - 4D-Thrive/Bliss | |
| # | - 4D-Trak | |
| # | - 2K Imaging | |
| # NMRF067 | SCANWISE IMPLANT | |
| Permanent | | 1303462/000020 |
| # NMRF090 | SMARTEXAM BRAIN | |
| Permanent | | 1303462/000020 |
| # NMRF112 | SMARTEXAM BREAST | |
| Permanent | | 1303462/000020 |
| # NMRF099 | SMARTEXAM KNEE | |
| Permanent | | 1303462/000020 |
| # NMRF098 | SMARTEXAM SHOULDER | |
| Permanent | | 1303462/000020 |
| # NMRF091 | SMARTEXAM SPINE | |
| Permanent | | 1303462/000020 |
| # NMRF077 | SPECTROSCOPY SPECIALIST | |
| Permanent | | 1303462/000020 |
| # NMRF118 | STARQUANT | |
| Permanent | | 1303462/000020 |
| # NMRF070 | SWI SPECIALIST | |
| Permanent | | 1303462/000020 |
| # NMRF123 | UPGR TO 4D-TRAK XD | |
| Permanent | | 1303462/000020 |
| # NMRF268 | Upgrade to latest sw version | |
| Permanent | | 1303462/000020 |

452213308911 IDL off-the-shelf software
31-Dec-9999 1303462/000020

989603040617 SW Key 17 K-T BLAST
31-Dec-9999 1303462/000020

989603040695 SW Key 95 mDIXON MRA
31-Dec-9999 1303462/000020

989603040734 SW Key 134 CS-SENSE Generic
31-Dec-9999 1303462/000020

989603040659 SW Key 59 Q-FLOW ANALYSIS
31-Dec-9999 1303462/000020

989603040652 SW Key 52 T1 PERFUSION
31-Dec-9999 1303462/000020

989603019801 LIC.MCAFEE VIRUSSCAN V8.5I EE
31-Dec-9999 1303462/000020

989603040636 SW Key 36 ADV PERIPHERAL ANGIO
31-Dec-9999 1303462/000020

989603040737 SW Key 137 CS-SENSE Body
31-Dec-9999 1303462/000020

989603040612 SW Key 12 ADVANCED EPI
31-Dec-9999 1303462/000020

989603040633 SW Key 33 MOVING BED IMA.
31-Dec-9999 1303462/000020

989603040705 SW Key 105 PATIENT CONDITIONS
31-Dec-9999 1303462/000020

989603040624 SW Key 24 PHASE NAVIGATORS
31-Dec-9999 1303462/000020

989603040682 SW Key 82 T2STAR R2STAR MAPPING
31-Dec-9999 1303462/000020

989603040614 Sw Key 14 3D BRAIN VIEW
31-Dec-9999 1303462/000020

989603040735 SW Key 135 CS-SENSE MSK
31-Dec-9999 1303462/000020

989603040643 SW Key 43 CARDIAC TRIGGERED IR
31-Dec-9999 1303462/000020

989603040710 SW Key 110 CENTRAKEYHOLE VIEW SHARING

31-Dec-9999 1303462/000020
989603040619 SW Key 19 EXTENDED FOV
31-Dec-9999 1303462/000020
989603040610 SW Key 10 3D BODY VIEW
31-Dec-9999 1303462/000020
989603019811 LIC.MACROMEDIA FLASH PLAYER10.1
31-Dec-9999 1303462/000020
989603040631 SW Key 31 CONTRAST ENH. PERF.
31-Dec-9999 1303462/000020
989603040605 SW Key 5 2D-VCG
31-Dec-9999 1303462/000020
989603019831 LIC.VXWORKS V5
31-Dec-9999 1303462/000020
989603040641 SW Key 41 MSMP
31-Dec-9999 1303462/000020
989603040645 SW Key 45 PROSP. MOTION CORR.
31-Dec-9999 1303462/000020
989603040691 SW Key 91 MULTI STATION VIEW
31-Dec-9999 1303462/000020
989603040692 SW Key 92 NEAR METAL VAT
31-Dec-9999 1303462/000020
989603040623 SW Key 23 DTI
31-Dec-9999 1303462/000020
989603040620 SW Key 20 SMART-SCAN SHOULDER
31-Dec-9999 1303462/000020
989603040637 SW Key 37 WHOLE BODY
31-Dec-9999 1303462/000020
989603040656 SW Key 56 DIFFUSION CALC.
31-Dec-9999 1303462/000020
989603040608 SW Key 8 2K IMAGING
31-Dec-9999 1303462/000020
989603040630 SW Key 30 CENTRAKEYHOLE
31-Dec-9999 1303462/000020
989603040647 SW Key 47 CARDIAC PERFUSION
31-Dec-9999 1303462/000020

989603040699 SW Key 99 AUTOVOICE
31-Dec-9999 1303462/000020

989603040634 SW Key 34 ADV. MOVING BED IMA.
31-Dec-9999 1303462/000020

989603040635 SW Key 35 ARTERIAL SPIN LABEL.
31-Dec-9999 1303462/000020

989603040640 SW Key 40 BODY RESP NAVIGATORS
31-Dec-9999 1303462/000020

989603040604 SW Key 4 SMART-SCAN KNEE
31-Dec-9999 1303462/000020

989603040708 SW Key 108 3D NERVEVIEW
31-Dec-9999 1303462/000020

989603040662 SW Key 62 MDIXON
31-Dec-9999 1303462/000020

989603040622 SW Key 22 DWI
31-Dec-9999 1303462/000020

989603040616 SW Key 16 HIGH SENSE FACTORS
31-Dec-9999 1303462/000020

989603040646 SW Key 46 REST GRID
31-Dec-9999 1303462/000020

989603040711 SW Key 111 3D BLACK BLOOD
31-Dec-9999 1303462/000020

989603040685 SW Key 85 XML/REC DATA EXPORT
31-Dec-9999 1303462/000020

989603040644 SW Key 44 RETROSPECTIVE
31-Dec-9999 1303462/000020

459800542671 Acronis ABR 11 OEM WS software license
31-Dec-9999 1303462/000020

989603040657 SW Key 57 PERFUSION CALC.
31-Dec-9999 1303462/000020

989603040606 SW Key 6 3D SPINE VIEW
31-Dec-9999 1303462/000020

989603040609 SW Key 9 3D MSK VIEW
31-Dec-9999 1303462/000020

989603040627 SW Key 27 FUNCTIONAL MRI

31-Dec-9999 1303462/000020
989603040611 SW Key 11 SPECTRO 3D
31-Dec-9999 1303462/000020
989603040693 SW Key 93 NEAR METAL SEMAC
31-Dec-9999 1303462/000020
989603019791 LIC.REAL VNC V4.2.6
31-Dec-9999 1303462/000020
989603040615 SW Key 15 SPECTRO SENSE
31-Dec-9999 1303462/000020
989603040650 SW Key 50 SPECTRO
31-Dec-9999 1303462/000020
989603040632 SW Key 32 CENTRAPLUSKEYHOLE
31-Dec-9999 1303462/000020
989603040719 SW Key 119 TSE OPTIMIZED REFOCUSING CONT
31-Dec-9999 1303462/000020
989603040648 SW Key 48 CARDIAC CAI
31-Dec-9999 1303462/000020
989603040625 SW Key 25 SMART-SCAN BREAST
31-Dec-9999 1303462/000020
989603040626 SW Key 26 TE>TR
31-Dec-9999 1303462/000020
989603040601 SW Key 1 SMART-SCAN SPINE
31-Dec-9999 1303462/000020
989603040677 SW Key 77 SWIP
31-Dec-9999 1303462/000020
989603040715 SW Key 115 INGENIA CONFIGURATION
31-Dec-9999 1303462/000020
989603040642 SW Key 42 MULTIVANE
31-Dec-9999 1303462/000020
989603040738 SW Key 138 CS-SENSE Cardiac
31-Dec-9999 1303462/000020
989603040621 SW Key 21 KEYHOLE
31-Dec-9999 1303462/000020
989603040638 SW Key 38 HIGH END RECON
31-Dec-9999 1303462/000020

#----- end of file -----

ECDb Version Information updated 11-May-2021 11:06:39 AM
Global ECDb Version : 00000000-00-00
Global ECDb installed Date : 11-MAY-2021 09:06:39
Regional ECDb : Not Installed

generation updated 7-May-2024 8:37:08 AM
Build Number 0083
Internal Release Number 99.99.99
Stream r571sp
Date 01-JUN-2023 19:28:04.78

Mrwin updated 11-May-2021 10:52:11 AM

Build Number 0172
Internal Release Number 99.99.99
Stream r571
Date 27-MAR-2020 14:19:11.75

MR R5 7 1 UserDocs14 updated 7-May-2024 8:37:17 AM
VERSION : 14
SUPPORTED_MR_RELEASE : 5.7.1
BUILD_NR. : V102.3.018
PRODUCT_MODEL : INGENIA
MANUFACTORY_VERSION : CEPNL
Installed at 07-May-2024 08:37:17

MR R5 7 1 UserDocs2 updated 11-May-2021 1:39:46 PM

VERSION : 2
SUPPORTED_MR_RELEASE : 5.7.1
BUILD_NR. : V102.2.011
PRODUCT_MODEL : INGENIA
MANUFACTORY_VERSION : CEPNL
Installed at 11-May-2021 13:39:46

R5 7 1 Servicepack1 updated 11-May-2021 10:40:25 AM

Release Number R5.7.1.1 (Build Number: 0023)
Installed at 11-MAY-2021 10:40:14

R5 7 1 Servicepack2 updated 9-Jun-2021 11:07:16 AM

Release Number R5.7.1.2 (Build Number: 0046)
Installed at 09-JUN-2021 11:07:05

R5 7 1 Servicepack3 updated 13-Feb-2022 7:38:11 AM

Release Number R5.7.1.3 (Build Number: 0057)
Installed at 13-FEB-2022 07:38:00

R5 7 1 Servicepack4 updated 7-May-2024 8:37:08 AM

Release Number R5.7.1.4 (Build Number: 0083)
Installed at 07-MAY-2024 08:36:57

Reconstructor Config updated 30-Jan-2026 4:07:29 AM

Platform: Windows 10 Enterprise 2016 LTSB 6.3, build:
14393.2125.amd64fre.rs1_release.180301-2139

Reconstructor clock speed: 2793 MHz
Reconstructor Model: HP Z420 Workstation
Reconstructor BIOS Version: J61 v03.52
Bulk Data Interface Driver: 6.1.7600.16385 built by: WinDDK Copyright (c) Philips
Electronics N.V. All rights reserved. INTEL BWin10-Recon version :
02.31.14.11 Date-time : 4:08 PM 3/28/201

Reconstructor0 has 4 processor(s) and 16384 Mbytes memory
Max memory available to reconstruction SW : 32768 Mbytes
Host computer model: HP Z420 Workstation
Host computer clock speed: 2800 Mhz
Host computer memory size: 34279931904 bytes
Host computer nr. of processors: 1
Host computer OS: Microsoft Windows 10 Enterprise 2016 LTSC
Host computer OS version: 10.0.14393
Host diskdrive model: StoreJet SCSI Disk Device, WD Elements 2621 USB
Device, WDC WD5000HHTZ-60N21V0, MTFDDAK128MAY-1AH1ZABHA

Driver Info

ACI 16 158037098001 451001024879 updated 3-Feb-2026 1:10:09 PM
Org. 12NC: 451001024879

Mseq: 0.
RXEser: 158036661034.
RXE12NC: 459800046225.

Mseq: 2.
RXEser: 158036661097.
RXE12NC: 459800046225.

Mseq: 1.
RXEser: 158024188068.

RXE12NC: 459800046225.

Mseq: 3.

RXEser: 158024188003.

RXE12NC: 459800046225.

ANTERIOR 00909 459801730112 updated 12-Jan-2026 3:01:10 PM
Org. 12NC: 459801730112

Mseq: 0.

RXEser: JS2408160077.

RXE12NC: 459800046226.

Mseq: 2.

RXEser: JS2408120456.

RXE12NC: 459800046226.

Mseq: 1.

RXEser: JS2408160086.

RXE12NC: 459800046226.

Mseq: 3.

RXEser: JS2408160081.

RXE12NC: 459800046226.

ANTERIOR 03559 459800118397 updated 3-Feb-2026 10:38:30 AM
Org. 12NC: 459800118397

Mseq: 0.

RXEser: JS2507180364.
RXE12NC: 459800046226.

Mseq: 2.
RXEser: JS2507180363.
RXE12NC: 459800046226.

Mseq: 1.
RXEser: JS2507180408.
RXE12NC: 459800046226.

Mseq: 3.
RXEser: JS2507180410.
RXE12NC: 459800046226.

CBH 171375-018 453567321856 updated 3-Feb-2026 1:10:08 PM
Org. 12NC: 453567321854

CDAS boards updated 30-Jan-2026 4:06:33 AM
Slot 1: empty
Slot 2: empty
Slot 3: empty
Slot 4: empty
Slot 5: empty
Slot 6: name: GRADM board id:1179587.3 fpga: v1 b1140
Slot 7: name: PAT board id:1179445.3 fpga: v15 b192
Slot 8: name: CBH board id:6732185.6 fpga: v3 b562
Slot 9: name: Clock board id:1179589.7 fpga: v17 b16
Slot 10: empty

Slot 11: empty
Slot 12: empty
Slot 13: name: Transmit1H board id:1179725.4 fpga: v1 b444
Slot 14: empty
Slot 15: empty
Slot 16: empty
Slot 17: empty
Slot 18: empty
Slot 19: empty

DCP 157027329007 459800443191 updated 28-Jan-2026 9:43:08 AM
Org. 12NC: 459800443191

Mseq: 0.
RXEser: 158127271019.
RXE12NC: 459800046226.

Mseq: 2.
RXEser: 158127271066.
RXE12NC: 459800046226.

Mseq: 1.
RXEser: 158127271074.
RXE12NC: 459800046226.

DCP RPB221773001 459800443191 updated 3-Feb-2026 1:10:09 PM
Org. 12NC: 459800443191

Mseq: 0.
RXEser: RPB185318385.

RXE12NC: 459800046225.

Mseq: 2.

RXEser: RPB185318305.

RXE12NC: 459800046225.

Mseq: 1.

RXEser: RPB185318370.

RXE12NC: 459800046225.

FLEX M ACI 3029 453530322322 updated 26-Jan-2026 6:40:59 AM
Org. 12NC: 453530322322

FLEX M ACI 347 459801729952 updated 26-Jan-2026 6:04:53 AM
Org. 12NC: 459801729952

FLEX S ACI 1803 452213264067 updated 2-Feb-2026 7:06:09 AM
Org. 12NC: 452213264067

FLEX S ACI 634 452213264067 updated 2-Feb-2026 6:31:28 AM
Org. 12NC: 452213264067

GRADIENT AMP updated 30-Jan-2026 4:06:36 AM
Rack 1

| | X | Y | Z | PS | ADCI |
|---------------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| NAME | 78X | 78X | 78X | 787 | ADCI |
| FIRMW VER | 01.50 | 01.50 | 01.50 | 00.88 | 00.66 |
| BOOT FIRMW VER | 00.98 | 00.98 | 00.98 | 00.98 | |
| SERIAL NO: FRU Ass. | 2119GS100. . 10 | 2314GS515. . 09 | 1116GS234. . 09 | | |
| | 1318GU020. . 12 | | | | |
| SERIAL NO: ModBoard | 151910380.0.E09 | 231416062.0.EF2 | 061613937.0.E07 | | |
| SERIAL NO: DrvBoard | 191914049.0.C04 | 492215833.0.C05 | 022212610.0.C05 | | |
| SERIAL NO: PowBoard | 191914386.0.C04 | 231416086.0.CC4 | 501519845.0.CC4 | | |
| SERIAL NO: FilBoard | 071315284.0.CC1 | 071315600.0.CC1 | 071315287.0.CC1 | | |
| SERIAL NO: IO Panel | 091316294.0.BD1 | 091316296.0.BD1 | 091316298.0.BD1 | | |
| SERIAL NO: OutFilB | | | | 121819022.0.B03 | |
| SERIAL NO: powDrvB | | | | 062426680.0.609 | |
| SERIAL NO: PSContB | | | | 071815692.0.303 | |
| SERIAL NO: 48 V PS | | | | 111818020.0.BB5 | |
| SERIAL NO: IO Board | | | | 091316456.b.1d2 | |
| SERIAL NO: GCI Bck | | | | | 061314909.0.BB3 |
| SERIAL NO: Rack | | | | | 1013GR01 .0.00H |
| VER INF EPLD_1 | 00.04 | 00.04 | 00.04 | 00.02 | |
| VER INF EPLD_2 | 00.02 | 00.02 | 00.02 | 00.08 | |
| VER INF EPLD_3 | 00.02 | 00.02 | 00.02 | | |
| VER INF DSP | 01.00 | 01.00 | 01.00 | | |

Rack 2

| | X | Y | Z | SwitchBox |
|---------------------|-----------------|-----------------|-----------------|-----------------|
| NAME | 78X | 78X | 78X | SBOX |
| FIRMW VER | 01.50 | 01.50 | 01.50 | 00.22 |
| BOOT FIRMW VER | 00.98 | 00.98 | 00.98 | 00.07 |
| SERIAL NO: FRU Ass. | 3911GS419. . 0E | 0609GS423. . 0B | 3817GS055. . 10 | |
| SERIAL NO: ModBoard | 371119456.0.EF2 | 060911263.0.608 | 281719209.0.E09 | |
| SERIAL NO: DrvBoard | 381119686.0.CC3 | 060911287.0.AC3 | 351712292.0.C04 | |
| SERIAL NO: PowBoard | 381119735.0.C0C | 060910882.0.B0B | 351712807.0.C04 | |
| SERIAL NO: FilBoard | 071315921.0.CC1 | 121819188.0.001 | 071315607.0.CC1 | |
| SERIAL NO: IO Panel | 091316281.0.BD1 | 091316295.0.BD1 | 091316301.0.BD1 | |
| SERIAL NO: ContBrd | | | | 041313780.0.CC2 |
| VER INF EPLD_1 | 00.04 | 00.04 | 00.04 | |

| | | | |
|----------------|-------|-------|-------|
| VER INF EPLD_2 | 00.02 | 00.02 | 00.02 |
| VER INF EPLD_3 | 00.02 | 00.02 | 00.02 |
| VER INF DSP | 01.00 | 01.00 | 01.00 |

GRADIENT GCI updated 30-Jan-2026 4:06:35 AM
IGCI

board model 4522 117 9719.7 fpga v3 b186

GRADIENT LCC updated 30-Jan-2026 4:06:35 AM

LCC version: 1.16

Serial: 10

Standby supported: yes

MAGNET MMU updated 3-Feb-2026 1:09:05 PM

MEU Version Number: 805:187.

Magnet type: WA15

NVC BASE 03493 459801290052 updated 3-Feb-2026 12:52:33 PM

Org. 12NC: 459801290052

Mseq: 0.

RXEser: 158152390048.

RXE12NC: 459800046226.

Mseq: 1.

RXEser: 158152390011.

RXE12NC: 459800046226.

NVC HEAD 05147 459800097954 updated 3-Feb-2026 12:38:15 PM
Org. 12NC: 459800097954

Mseq: 0.
RXEser: 158010622032.
RXE12NC: 459800046225.

Mseq: 1.
RXEser: 158010622092.
RXE12NC: 459800046225.

NVC HEAD NECK 06436 459800097965 updated 3-Feb-2026 3:09:23 AM
Org. 12NC: 459800097965

Mseq: 0.
RXEser: 158068989088.
RXE12NC: 459800046225.

Mseq: 1.
RXEser: 158068989030.
RXE12NC: 459800046225.

PATCOM - AudioModule updated 30-Jan-2026 4:04:47 AM
AudioModule 12 NC : 459801497442
AudioModule SW Version : 0x0128
AudioModule SW Build : 0x1307
AudioModule HW Version : AM3_AM5_V20

AudioModule Serial number : 342057

PATCOM - Host-AIBo updated 30-Jan-2026 4:04:45 AM
Host-AIBo 12 NC : 459800138513
Host-AIBo SW Version : 0x0128
Host-AIBo SW Build : 0x1307
Host-AIBo HW Version : AIBo_020_v01
Host-AIBo Serial number : 006396

PATSUP-MT PLC updated 30-Jan-2026 4:04:49 AM
PLC extension 12NC: 45221179729.9
PLC extension ID: 8
PLC extension SW version: 9
PLC firmware version: 2
PLC firmware build: 12
PLC TSW-part: 11
PLC hardware type: 3

PHYSIOLOGY PPU Sensor version(last used) updated 31-Jan-2026 12:05:58 AM
wPPU version = 0.23

PHYSIOLOGY updated 30-Jan-2026 4:06:33 AM
wBTU version: FW = 10.2; SW = 3.2; wECG = N/A; wPPU = N/A

POSTERIOR 00035 452213277806 updated 3-Feb-2026 1:10:09 PM
Org. 12NC: 452213277806

PSU IRF2 Device 459800505121 updated 3-Feb-2026 1:10:09 PM
Org. 12NC: 459800505121

RF-TX 1H RF-AMP updated 30-Jan-2026 4:06:35 AM

Amp. Channel 1 (CDAS.Slot13)
RF amplifier model: 53-S30-64
RF amplifier serial: PS09381581
RF amplifier total hours: 42544
RF amplifier version info:

MKS Instruments
63.87 MHz Solid State Amplifier

SW P/N SRO-X336 Rev A CPLD Ver 0
SW P/N SRO-X332 Rev I IOP Ver 20/01/2012 Hours: 11
SW P/N SRO-X337 Rev B IOPBL Ver 23/03/2007 Hours: 11
SW P/N SRO-X331 Rev H RFP Ver 02/09/2014 Hours: 10
SW P/N SRO-X340 Rev B PSC Ver 06/01/2006 Hours: 07

SBM Device updated 1-Feb-2026 2:59:05 PM
SBM firmware version: 5
SBM firmware build : 2058

SBM S181500299 452211800044 updated 3-Feb-2026 1:10:08 PM
Org. 12NC: 452211800044

SENSE ANKLE 8 ACI 000370 453530252261 updated 1-Feb-2026 2:53:48 AM
Org. 12NC: 453530252261

SENSE BREAST 16 ACI 001131 453530280721 updated 3-Feb-2026 5:59:58 AM
Org. 12NC: 453530280721

SENSE KNEE 16 ACI 000590 453530252671 updated 2-Feb-2026 12:07:21 PM
Org. 12NC: 453530252671

SENSE SHOULDER8 DCI 000155 453530283042 updated 2-Feb-2026 3:03:03
PM
Org. 12NC: 453530283042

Mseq: 0.

RXEser: PPB504927296.

RXE12NC: 459800046225.

Mseq: 1.

RXEser: PPB504927292.

RXE12NC: 459800046225.

SENSE SMALL EXTR 8 ACI 000312 453530283161 updated 2-Feb-2026

10:47:26 PM

Org. 12NC: 453530283161

SENSE WRIST 8 ACI 000225 453530252281 updated 31-Jan-2026 8:29:36 PM

Org. 12NC: 453530252281

TFINT PPB504194007 452213700426 updated 3-Feb-2026 1:10:09 PM

Org. 12NC: 452213700426

Mseq: 0.

RXEser: PPB504158303.

RXE12NC: 459800046225.

XRUI updated 30-Jan-2026 4:04:47 AM

System Control Unit-3 ,front left 12NC nr(last 8 dig) 1331397.4 ,version 2 ,build 22

System Control Unit-3 ,front right 12NC nr(last 8 dig) 1331396.4 ,version 2 ,build 22

User Interface Hub 12NC nr(last 8 dig) 1179731.5 ,boot version 1 ,build 16

UIH APPL software version 1 ,build 61

UIH FPGA software version 0 ,build 192

UIH CPLD software version 0 ,build 128

Install Information

Philips Applications

Version

Install Date

Philips ICE Client x64 3.3.0.20190703 3.3.0.20190703 2024-05-07 19:21

Philips Healthcare Fussite Client Package-3.5.279.1534 3.5.279.1534 2021-05-12
13:38

Other Applications

Version

Install Date

Windows Driver Package - X-Rite (Chroma) XRiteDevices (08/21/2006 2.0.0.0)
08/21/2006 2.0.0.0 2021-05-11 11:07

Service Pack 2 for SQL Server 2016 (KB4052908) (64-bit) 13.2.5026.0 2021-05-11
10:37

Microsoft SQL Server 2016 (64-bit) 2019-12-02 10:52

Microsoft SQL Server 2016 (64-bit) 2019-12-02 10:52

Windows Media Encoder 9 Series x64 Edition 2019-12-02 10:56

SQL Server 2016 Database Engine Services 13.2.5026.0 2021-05-12 13:38

Sql Server Customer Experience Improvement Program 13.2.5026.0 2021-05-12
13:38

Microsoft SQL Server 2016 Setup (English) 13.2.5026.0 2021-05-12 13:38

Microsoft SQL Server 2012 Native Client 11.3.6518.0 2021-05-12 13:38

SQL Server 2016 Common Files 13.2.5026.0 2021-05-12 13:38

VNC Server 6.1.0 6.1.0.27437 2021-05-12 13:38

PMS MR Interfaces 33.24.172.20200327152215 mrmmain 2021-05-12 13:38

PMS MR Platform 33.24.172.20200327152229 mrmmain 2024-05-07 19:21

Microsoft Visual C++ 2010 x64 Redistributable - 10.0.40219 10.0.40219
2021-05-12 13:38

7-Zip 18.05 (x64 edition) 18.05.00.0 2021-05-12 13:38

MR_R5_7_1_ServicePack3 6.1.514.20211004060113 mrmmain 2024-05-07 19:21

SQL Server 2016 DMF 13.0.1601.5 2021-05-12 13:38

SQL Server 2016 DMF 13.0.1601.5 2021-05-12 13:38

PMS MR FieldService 33.24.172.20200327152208 mrmmain 2024-05-07 19:21

McAfee Agent 5.5.0.447 2024-05-07 19:21

Microsoft VSS Writer for SQL Server 2016 13.2.5026.0 2021-05-12 13:38

SQL Server 2016 Shared Management Objects 13.0.16107.4 2021-05-12 13:38

PMS MR Scansoftware 33.24.172.20200327152125 mrmmain 2024-05-07 19:21

Microsoft SQL Server 2016 T-SQL ScriptDom 13.2.5026.0 2021-05-12 13:38

SQL Server 2016 Connection Info 13.0.16108.4 2021-05-12 13:38

SQL Server 2016 Database Engine Services 13.2.5026.0 2024-05-07 19:21

SQL Server 2016 Common Files 13.2.5026.0 2024-05-07 19:21

Microsoft Visual C++ 2008 Redistributable - x64 9.0.30729.6161 9.0.30729.6161
2021-05-12 13:38

Microsoft SQL Server 2008 Setup Support Files 10.3.5500.0 2024-05-07 19:21

MR_R5_7_1_ServicePack1 6.1.514.20200721103856 mrmmain 2021-05-12 13:38

SQL Server 2016 Database Engine Shared 13.2.5026.0 2024-05-07 19:21

PMS MR DeviceControl 33.24.172.20200327151858 mrmmain 2024-05-07 19:21

SQL Server 2016 Connection Info 13.0.16108.4 2021-05-12 13:38

SQL Server 2016 SQL Diagnostics 13.0.1601.5 2021-05-12 13:38

MR_R5_7_1_ServicePack2 6.1.514.20210302051834 mrmmain 2024-05-07 19:21

Microsoft SQL Server 2016 RsFx Driver 13.2.5026.0 2021-05-12 13:38

SQL Server 2016 Database Engine Shared 13.2.5026.0 2021-05-12 13:38

SQL Server 2016 XEvent 13.0.1601.5 2021-05-12 13:38

PSC 33.24.172.20200327151040 mrmmain 2021-05-12 13:38

IDL 7.0.6 7.0.6.0000 2021-05-12 13:38

Microsoft ODBC Driver 13 for SQL Server 13.2.5026.0 2021-05-12 13:38

Microsoft Visual C++ 2005 Redistributable (x64) 8.0.61000 2021-05-12 13:38

MR_R5_7_1_ServicePack4 6.1.514.20230601201449 mrmmain 2024-05-07 19:21

SQL Server 2016 Shared Management Objects 13.0.16107.4 2021-05-12 13:38

VNC Viewer 6.1.0 6.1.0.27437 2021-05-12 13:38

NVIDIA Ansel 7.0.548.0 2021-05-12 13:38

NVIDIA Control Panel 441.28 441.28 2021-05-12 13:38

NVIDIA Graphics Driver 441.28 441.28 2021-05-12 13:38

NVIDIA nView 149.77 149.77 2021-05-12 13:38

NVIDIA WMI 2.34.0 2.34.0 2021-05-12 13:38

NVIDIA HD Audio Driver 1.3.38.21 1.3.38.21 2021-05-12 13:38

NVIDIA Install Application 2.1002.334.1 2021-05-12 13:38

NVIDIA Display Container 1.19 2021-05-12 13:38

NVIDIA Display Container LS 1.19 2021-05-12 13:38

NVIDIA Display Watchdog Plugin 1.19 2021-05-12 13:38

NVIDIA Display Session Container 1.19 2021-05-12 13:38

SQL Server 2016 Shared Management Objects Extensions 13.2.5026.0
2021-05-12 13:38

PMS MR Protocols 33.24.172.20200327131318 mrmmain 2024-05-07 19:21

PMS MR ImageHandling 33.24.172.20200327152246 mrmmain 2024-05-07 19:21

PMS MR DAS 33.24.172.20200327145029 mrmmain 2021-05-12 13:38

Windows Media Encoder 9 Series x64 Edition 10.0.0.3809 2021-05-12 13:38

Mrwin 33.24.172.20200327143536 mrmmain 2021-05-12 13:38

SQL Server 2016 Batch Parser 13.0.1601.5 2021-05-12 13:38

PMS MR Coildata 33.24.172.20200327152228 mrmmain 2021-05-12 13:38

PMS MR Recon 33.24.172.20200327152108 mrmmain 2024-05-07 19:21

SQL Server 2016 XEvent 13.0.1601.5 2021-05-12 13:38

Microsoft Visual C++ 2015 x64 Additional Runtime - 14.0.24212 14.0.24212
2021-05-12 13:38

SQL Server 2016 Full text search 13.2.5026.0 2024-05-07 19:21

SQL Server 2016 Shared Management Objects Extensions 13.2.5026.0
2021-05-12 13:38

Microsoft Visual C++ 2015 x64 Minimum Runtime - 14.0.24212 14.0.24212
2021-05-12 13:38

MediCal QAWeb Agent 1.13.1400 2024-05-07 19:21

Microsoft SQL Server 2016 T-SQL Language Service 13.0.14500.10 2021-05-12
13:38

Copssh 5.9.0 (remove only) 2021-05-11 10:49

IDL 7.0.6 2021-05-11 10:37

Strawberry Perl 5.14.2001 2021-05-12 13:38

Acronis Backup 11.7 Agent for Windows 11.7.44421 2024-05-07 19:21

Acronis Backup 11.7 Command-Line Tool 11.7.44421 2021-05-12 13:38

Microsoft Visual C++ 2015 Redistributable (x64) - 14.0.24212 14.0.24212.0
2021-05-11 10:37

Microsoft Visual C++ 2015 x86 Minimum Runtime - 14.0.24212 14.0.24212
2021-05-12 13:38

Microsoft Visual C++ 2015 Redistributable (x86) - 14.0.24212 14.0.24212.0
2021-05-11 10:37

Microsoft SQL Server Data-Tier Application Framework (x86) 13.0.3225.4
2021-05-12 13:38

Browser for SQL Server 2016 13.2.5026.0 2021-05-12 13:38

Microsoft Visual C++ 2005 Redistributable 8.0.61001 2021-05-12 13:38

Microsoft Visual C++ 2015 x86 Additional Runtime - 14.0.24212 14.0.24212
2021-05-12 13:38

Microsoft Visual C++ 2008 Redistributable - x86 9.0.30729.6161 9.0.30729.6161
2021-05-12 13:38

Adobe Acrobat Reader DC 19.021.20049 2021-05-12 13:38

McAfee VirusScan Enterprise 8.8.013000 2021-05-12 13:38

Microsoft Visual C++ 2010 x86 Redistributable - 10.0.40219 10.0.40219
2021-05-12 13:38

Realtek High Definition Audio Driver 6.0.1.7161 2021-05-12 13:38

Acronis Backup 11.7 Agent Core 11.7.44421 2024-05-07 19:21

Version check of Microsoft Windows

7 all 7 files checked OK

Version check of PMS MR Coildata

904 all 904 files checked OK

Version check of PMS MR DAS

9 all 9 files checked OK

Version check of PMS MR DeviceControl

862 all 862 files checked OK

Version check of PMS MR FieldService

7630 all 7630 files checked OK

Version check of PMS MR ImageHandling

7443 all 7443 files checked OK

Version check of PMS MR Platform

2625 all 2625 files checked OK

Version check of PMS MR Protocols

15 all 15 files checked OK

Version check of PMS MR Recon

618 all 618 files checked OK

Version check of PMS MR Scansoftware

3866 all 3866 files checked OK
Version check of PSC

70 all 70 files checked OK
Version check of Debugger / Development Tools
Development tools should not be present on a system.

5 all 5 files checked OK

SCSI Devices

Scsi Port 2

Interrupt = 0x0000000B (11)
Dma64BitAddresses = 0x00000001 (1)
Driver = iaStorA

Scsi Bus 0

Initiator Id 255

Target Id 0

Logical Unit Id 0

InquiryData REG_?????(3)
SerialNumber = 14310CEE86AC
Identifier = ATA MTFDDAK128MAY-1AM504
DeviceType = DiskPeripheral
DeviceIdentifierPage REG_?????(3)

Scsi Bus 1

Initiator Id 255

Target Id 0

Logical Unit Id 0

InquiryData REG_?????(3)
SerialNumber = WD-WXT1E64JL38A
Identifier = ATA WDC WD5000HHTZ-66A02
DeviceType = DiskPeripheral
DeviceIdentifierPage REG_?????(3)

Scsi Bus 2

Initiator Id 255

Target Id 0

Logical Unit Id 0

InquiryData REG_?????(3)
Identifier = hp DVD A DH16AESH 3HH5
DeviceType = CdRomPeripheral
Scsi Bus 3
Initiator Id 255
Scsi Bus 4
Initiator Id 255
Scsi Bus 5
Initiator Id 255
Scsi Bus 6
Initiator Id 255
Scsi Bus 7
Initiator Id 255
Scsi Bus 8
Initiator Id 255
Scsi Port 3
Dma64BitAddresses = 0x00000001 (1)
Driver = UASPStor
Scsi Bus 0
Initiator Id 255
Target Id 0
Logical Unit Id 0
InquiryData REG_?????(3)
SerialNumber = 807ZA0CE29XW-DW
Identifier = StoreJet 0
DeviceType = DiskPeripheral
DeviceIdentifierPage REG_?????(3)

Microsoft Diagnostics Report

Microsoft Windows platform: NT/2000/XP

Version 6.2 build 9200

• System Information report written at: 02/03/26 13:15:50

System Name: PHILIPS-QKTV6EF

[System Summary]

Item Value

OS Name Microsoft Windows 10 Enterprise 2016 LTSC

Version 10.0.14393 Build 14393

Other OS Description Not Available

OS Manufacturer Microsoft Corporation

System Name PHILIPS-QKTV6EF

System Manufacturer Hewlett-Packard

System Model HP Z420 Workstation

System Type x64-based PC

System SKU LJ449AV

Processor Intel(R) Xeon(R) CPU E5-1603 0 @ 2.80GHz, 2800 Mhz, 4 Core(s), 4 Logical Processor(s)

BIOS Version/Date Hewlett-Packard J61 v03.88, 5/27/2015

SMBIOS Version 2.7

Embedded Controller Version 255.255

BIOS Mode Legacy

BaseBoard Manufacturer Hewlett-Packard

BaseBoard Model Not Available

BaseBoard Name Base Board

Platform Role Workstation

Secure Boot State Unsupported

PCR7 Configuration Binding Not Possible

Windows Directory C:\Windows

System Directory C:\Windows\system32

Boot Device \Device\HarddiskVolume1

Locale United States

Hardware Abstraction Layer Version = "10.0.14393.3297"

User Name PHILIPS-QKTV6EF\MRService

Time Zone Jerusalem Standard Time

Installed Physical Memory (RAM) 32.0 GB

Total Physical Memory 31.9 GB

Available Physical Memory 20.9 GB

Total Virtual Memory 47.9 GB

Available Virtual Memory 33.9 GB

Page File Space 16.0 GB

Page File G:\pagefile.sys

Device Guard Virtualization based security Not enabled

Hyper-V - VM Monitor Mode Extensions Yes

Hyper-V - Second Level Address Translation Extensions Yes

Hyper-V - Virtualization Enabled in Firmware No

Hyper-V - Data Execution Protection Yes

[Hardware Resources]